



## **Analytical Balance TRAB-612**

[www.safetyzon.com](http://www.safetyzon.com) | [info@safetyzon.com](mailto:info@safetyzon.com)

## Description

Analytical Balance TRAB-612 is a fully automatic, internal calibration device with a completely transparent glass windshield. It has a huge backlit LCD, which improves the operator's visibility. It can weigh up to 210 g from 0 and has a minimum reading of 0.001 g. It is a faster and much more reliable device with a double CAM experiment.

## Features:

- > It has humanized design, novel structure, and 5 fully transparent glass windshield surfaces
- > With the piece, unit conversion, full-range peel, zero tracking, and other features such as multi-functionality, simple and dependable operation
- > Includes an RS232C output interface that can be directly connected to a computer, printer, another device
- > It employs a double CAM experiment with a fully automatic internal calibration function, which increases its stability and speed
- > The Balance includes applications such as preheating, calibration, screen saver, and others
- > The large backlit LCD screen improves the operator's vision by making it more soft and clear

## Specifications :

<b>Linear Error</b>	±0.0005 g
<b>Minimum Reading</b>	0.0001 g
<b>Output Interface</b>	Rs232C
<b>Pan Size</b>	Ø80 mm
<b>Power</b>	220V+22V/-33V, 50Hz
<b>Repeatability Error</b>	0.0002 g
<b>Scale Interval</b>	0.1 mg
<b>Size</b>	330×215×350 mm
<b>Weighing Range</b>	0 to 210 g
<b>Weight</b>	7.5 Kg



West 110th Street Overland Park KS 66210, USA  
Email: [info@safetyzon.com](mailto:info@safetyzon.com) | Website: [www.safetyzon.com](http://www.safetyzon.com)