



## **Automatic Autoclave TRA-601**

[www.safetyzon.com](http://www.safetyzon.com) | [info@safetyzon.com](mailto:info@safetyzon.com)

## Description

Automatic Autoclave TRA-601 is a microprocessor-controlled tabletop electronic device with fully automatic, pre-programmed sterilization cycles and an LCD screen that displays continuous cycle information. It has an 8 L capacity, a maximum temperature of 134 °C, and a pressure of 210 MPa. It has an electro-polished stainless steel chamber and door. A simple user interface is very easy to use.

## Features:

- > Tabletop autoclaves that are full auto and have pre-programmed sterilization cycles
- > The chamber and entrance are made of stainless steel and have an electro-polish finish
- > Control system with high precision for perfect sterilization results
- > The built-in independent steam generator ensures that there is always steam accessible for faster cycles
- > The user interface is very friendly and can operate very easily
- > The LCD screen constantly provides data about the cycle
- > The open-type water tank is very simple to fill and clean
- > Microprocessor-controlled electronic autoclave that is fully automated
- > It is a thermal autoclave with a non-vacuum working method
- > Each operating status and program could be displayed as an icon on the LCD screen

## Specifications :

<b>Capacity</b>	8 L
<b>Dimension</b>	370×345×565 mm
<b>Heating Power</b>	1400 W
<b>Maximum Pressure</b>	210 MPa
<b>Maximum Temperature</b>	134 °
<b>Maximum Timer</b>	0 to 55 min
<b>Number Of Trays</b>	2
<b>Power Supply</b>	220V 50Hz /110V60Hz
<b>Sterilization Room Size</b>	Ø170×320
<b>Weight</b>	35 kg



West 110th Street Overland Park KS 66210, USA  
Email: [info@safetyzon.com](mailto:info@safetyzon.com) | Website: [www.safetyzon.com](http://www.safetyzon.com)