



## Automatic Sterilizer TRSA-602

[www.safetyzon.com](http://www.safetyzon.com) | [info@safetyzon.com](mailto:info@safetyzon.com)

## Description

Automatic Sterilizer TRSA-602 is a Computer Control Auto Recycle Sterilizer with a Stainless Steel Design. It has a capacity of 50 L and a maximum pressure range of 0.145 to 0.165Mpa. It delivers accurate results within a temperature range of 50°C to 126°C. Equipped with a quick-opening door framework with a hand wheel and a mechanical safety-interlocking device.

## Features:

- > It has a quick-open door structure with a hand wheel and a mechanical safety-interlocking device
- > A computer-controlled auto-recycled sterilization mechanism is built into the device
- > The LED screen displays the current state of operation
- > The temperature range is 50°C to 126°C
- > Having with the period of 0 to 9999min
- > Auto-discharging due to overpressure 0.145 to 0.165 mmHg
- > Sterilization is automatically terminated with a beep
- > Stainless steel is the only material used in its construction

## Specifications :

<b>Capacity</b>	50 L
<b>Dimension</b>	51×51×125cm
<b>Heating Power</b>	3.5 KW
<b>Maximum Pressure</b>	0.145-0.165Mpa
<b>Maximum Temperature</b>	126 ?
<b>Maximum Timer</b>	9999min
<b>Power Supply</b>	220V 50Hz
<b>Sterilization Room Size</b>	?340×425mm
<b>Weight</b>	58 kg



West 110th Street Overland Park KS 66210, USA  
Email: [info@safetyzon.com](mailto:info@safetyzon.com) | Website: [www.safetyzon.com](http://www.safetyzon.com)