



Automatic Sterilizer TRSA-603

www.safetyzon.com | info@safetyzon.com

Description

Automatic Sterilizer TRSA-603 is a Computer Control Auto Recycle Sterilizer with a Stainless Steel Design. It has a capacity of 75 L and a maximum pressure range of 0.145 to 0.165Mpa. It delivers accurate results within a temperature range of 50°C to 126°C. Equipped with a quick-opening door framework with a hand wheel and a mechanical safety-interlocking device.

Features:

- > It has a quick-open door structure with a hand wheel and a mechanical safety-interlocking device
- > A computer-controlled auto-recycled sterilization mechanism is built into the device
- > The LED screen displays the current state of operation
- > The temperature range is 50°C to 126°C
- > Having with the period of 0 to 9999min
- > Auto-discharging due to overpressure 0.145 to 0.165 mmHg
- > Sterilization is automatically terminated with a beep
- > Stainless steel is the only material used in its construction

Specifications :

Capacity	75 L
Dimension	57×57×130 cm
Heating Power	4.5 KW
Maximum Pressure	0.145-0.165Mpa
Maximum Temperature	126 ?
Maximum Timer	9999min
Power Supply	380V 50Hz
Sterilization Room Size	?340×425mm
Weight	77 kg



West 110th Street Overland Park KS 66210, USA
Email: info@safetyzon.com | Website: www.safetyzon.com