



## Balance TRB-602

[www.safetyzon.com](http://www.safetyzon.com) | [info@safetyzon.com](mailto:info@safetyzon.com)

## **Description**

Balance TRB-602 is a compact device with high stability; the sensor uses a unique locking device with ac/dc level, an LCD backlight display, and a weighing reaction speed range of 0 to 1000 g with precision of 10 mg. It is intended to be used in conjunction with a printer or computer.

## **Features:**

- > It has ac/dc, level, LCD backlight display, fast weighing reaction time, and high stability
- > The sensor employs a one-of-a-kind locking mechanism
- > Counting function, percentage, automatic calibration, multiple calibration options, calibration weights
- > With unit conversion, stability, memory function, animal weighing, full range gross weight reduction, and so on
- > With an RS232 interface for connecting to a printer or computer

## Specifications :

|                        |                |
|------------------------|----------------|
| <b>Accuracy</b>        | 10 mg          |
| <b>Minimum Reading</b> | 0.01g          |
| <b>Pan Size</b>        | ∅130 mm        |
| <b>Range</b>           | 0 to 1000 g    |
| <b>Size</b>            | 310×235×150 mm |
| <b>Weight</b>          | 1 Kg           |



West 110th Street Overland Park KS 66210, USA  
Email: [info@safetyzon.com](mailto:info@safetyzon.com) | Website: [www.safetyzon.com](http://www.safetyzon.com)