



## Balance TRB-603

[www.safetyzon.com](http://www.safetyzon.com) | [info@safetyzon.com](mailto:info@safetyzon.com)

## **Description**

Balance TRB-603 is a compact device with high stability; the sensor uses a unique locking device with ac/dc level, an LCD backlight display, and a weighing reaction speed range of 0 to 2000 g with a precision of 10 mg. It is intended to be used in conjunction with a printer or computer.

## **Features:**

- > It has ac/dc, level, LCD backlight display, fast weighing reaction time, and high stability
- > The sensor employs a one-of-a-kind locking mechanism
- > Counting function, percentage, automatic calibration, multiple calibration options, calibration weights
- > With unit conversion, stability, memory function, animal weighing, full range gross weight reduction, and so on
- > With an RS232 interface for connecting to a printer or computer

## Specifications :

<b>Accuracy</b>	10 mg
<b>Minimum Reading</b>	0.01g
<b>Pan Size</b>	∅130 mm
<b>Range</b>	0 to 2000 g
<b>Size</b>	310×235×150 mm
<b>Weight</b>	1 Kg



West 110th Street Overland Park KS 66210, USA  
Email: [info@safetyzon.com](mailto:info@safetyzon.com) | Website: [www.safetyzon.com](http://www.safetyzon.com)