



Digital Heating Water Bath TRHB-601

www.safetyzon.com | info@safetyzon.com

Description

Digital Heating Water Bath TRHB-601 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °. It holds with a capacity of 3 L. It is designed with a stainless steel interior with an aluminum ring cover.

Features:

- > It is equipped with a microprocessor control
- > mechanism with a stainless steel interior
- > Designed with a digital display inner temperature
- > It is covered with an aluminum ring
- > It is a thermostatic device

Specifications :

Capacity	3 L
Dimension	23x23x24 cm
Power	300 W
Temperature Fluctuation	± 0.5 ?
Temperature Homogeneity	± 1 ?
Temperature Range	(RT+5) -300 ?
Voltage	110V 60Hz /220V 50 Hz
Weight	3.5kg
Working Size	16x16x12 cm



West 110th Street Overland Park KS 66210, USA
Email: info@safetyzon.com | Website: www.safetyzon.com