



Digital Heating Water Bath TRHB-603

www.safetyzon.com | info@safetyzon.com

Description

Digital Heating Water Bath TRHB-603 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °. It holds with a capacity of 10 L. It is designed with a stainless steel interior with an aluminum ring cover.

Features:

- > It is equipped with a microprocessor control
- > mechanism with a stainless steel interior
- > Designed with a digital display inner temperature
- > It is covered with an aluminum ring
- > It is a thermostatic device

Specifications :

Capacity	10 L
Dimension	37x37x24 cm
Power	1000 W
Temperature Fluctuation	±1 °
Temperature Homogeneity	±1 °
Temperature Range	(RT+5) -300 °
Voltage	110V 60Hz /220V 50 Hz
Weight	5.5 kg
Working Size	30x30x12 cm



West 110th Street Overland Park KS 66210, USA
Email: info@safetyzon.com | Website: www.safetyzon.com