



## Multipurpose Horizontal Electrophoresis Cell TRME-601

[www.safetyzon.com](http://www.safetyzon.com) | [info@safetyzon.com](mailto:info@safetyzon.com)

## Description

Multipurpose Horizontal Electrophoresis Cell TRME-601 is primarily used for agarose gels of DNA and RNA, and it can add samples via pipette. It includes the Multi-specification gel tray: 30×130mm; 130×65mm; 65×130mm; 65×65mm, as well as PCR electrophoresis with the comb, 27 wells/4 rows. It has a moveable e-engine and is designed without a washer. It is a straightforward replacement with no leaks.

## Features:

- > PCR electrophoresis: 27 wells\*4 rows, 108 samples obtainable (including maker)
- > Built-in movable e-engine with no washer, simple replacement, and no leakage
- > Gel tray with multiple specifications: 30130mm; 13065mm; 65130mm; 6565mm
- > The comb with 14, 19, and 27 wells supports 8 and 12 sample adding channels, respectively
- > It can add through a pipette and is primarily used for agarose gels of DNA and RNA

## Specifications :

<b>Buffering liquid volume</b>	1000ml
<b>Comb specification</b>	0.75mm:7+7 wells/14 wells, 9+9 wells/19wells, 1.0 mm: 2.0mm: 3+2 wells/3+3 wells12+12 wells/27 wells, 1.5mm: 7+7 wells/14 wells, 9+9 wells/19 wells;
<b>Dimension (L×W×H)</b>	300×170×80mm
<b>Gel area</b>	130×130mm, 130×65mm, 65×130mm, 65×65mm
<b>Weight</b>	1.1Kg



West 110th Street Overland Park KS 66210, USA  
Email: [info@safetyzon.com](mailto:info@safetyzon.com) | Website: [www.safetyzon.com](http://www.safetyzon.com)