



NC Linear Shaker TRLS-601

www.safetyzon.com | info@safetyzon.com

Description

NC Linear Shaker TRLS-601 has dual LCD screens for parameter setting and viewing. It can withstand a maximum load weight and operates in a timed and constantly functioning mode. The dual LCD screens separately display time and speed, with a speed range of 40 to 200 rpm and a time range of 1 minute to 19 hours 59 minutes.

Features:

- > There is a maximum load-bearing weight of 7.5kg allowed
- > The period is from 1min 19 hrs. 59 min
- > The time and speed are displayed separately on dual LCD screens
- > There are numerous inserts and fittings available
- > Designed with a brushless DC motor requires no maintenance

Specifications :

Ambient Temperature	5 to 40°C
Amplitude	20mm
Humidity	0.8
Maximum Load	3kg
Motor Type	DC motor
Speed Display	LCD
Speed Range	40 to 200rpm
Time setting range	0 to 19h 59min
Voltage	100 to 240V, 50/60Hz



West 110th Street Overland Park KS 66210, USA
Email: info@safetyzon.com | Website: www.safetyzon.com