



Pharmaceutical Refrigerator TRPR-601

www.safetyzon.com | info@safetyzon.com

Description

Pharmaceutical Refrigerator TRPR-601 has a microprocessor-controlled digital display and temperature adjustment with a 0.1°C increment. For simple inspection, a huge digital display monitor is used. Presented with a cutting-edge forced cooling system to keep the interior temperature stable and uniform it has a comprehensive temperature alarm system. Temperature variation of less than 3 °C with an ambient temperature of 10 to 32 °C. Lock to prevent unauthorized access, forced-air cooling with an optimized air distribution system.

Features:

- > Special compressor with high efficiency and field-proven reliability A permanently lubricated cooling fan for safety and longevity
- Forced-air cooling with an optimized air distribution system is designed To attain maximum temperature consistency and stability
- Optimized refrigeration system intention for more actual cooling and quick retrieval Temperature variations are limited to 3 °C A temperature range of 2 to 8 °C is adjustable Temperature adjustment in 0.1 °C increments and a digital display controlled by a microprocessor large digital display screen Safety door lock intention to avoid door open Big digital screen display for easy observation Optional testing hole Adjustable shelves Inner light design for a clear display
- > Designed with a complete temperature alarm system featuring an audible buzzer and a visible flashing light
- > Failures due to high and low-temperature sensor error, door opening, power failure, and the low battery will give an alert
- > Progressed constrained cooling system to keep within temperature steady and uniform
- > Evaporator with an uncommon cycle plan of the constrained air cooling system to guarantee no-ice inside
- > Advanced space use plan with staggered and customizable rack tallness to oblige capacity necessity of drug results of different sorts
- > Wellbeing lock to forestall unapproved access
- > Caster plan and even out legs LED Interior light
- > Having with auto defrost
- > The material is color sprayed with a PURF insulation method

Specifications :

Ambient Temperature	10 to 32 ?
Capacity	60/2.11(L/cu. Ft)
Current	0.59 A
Exterior Size(W×D×H)	495×450×710 mm
G. Wt.	37 kg
Interior Size(W×D×H)	410×315×560 mm
N. Wt.	33 kg
Package Size(W×D×H)	558×510×760 mm
Power	72 W
Refrigerant	R600a
Shelf / Drawer	Shelves / 2
Temperature Range	2 to 8 ?
Voltage/Frequency	220/50 Hz



West 110th Street Overland Park KS 66210, USA
Email: info@safetyzon.com | Website: www.safetyzon.com