



Portable Sterilizer TRPS-604

www.safetyzon.com | info@safetyzon.com

Description

Portable Sterilizer TRPS-604 is a stainless steel device with an electric/coal heating mode and a Dual scale numerical indication pressure gauge. It is having with a movable interior that is convenient for operation and with a capacity of 24 L and disinfecting Pressure of 0.14-0.165Mpa. It is having with the self-inflating type seal.

Features:

- > It is made of stainless steel
- > It has an overpressure auto-discharging system. 0.145 to 0.165MPa
- > The maximum working temperature is 126 °C
- > Numeric indication with a dual magnitude Seal that self-inflates
- > There are two modes: electric heating and electric/coal heating
- > The movable interior makes operation easier

Specifications :

Capacity	24 L
Dimension	39×39×54 cm
Disinfecting Pressure	0.14-0.165 MPa
Heating Mode	Electric/coal heating
Maximum Working Temperature	126 ?
Power	2 KW
Voltage	110 V 60 Hz /220 V 50 Hz
Weight	17 kg



West 110th Street Overland Park KS 66210, USA
Email: info@safetyzon.com | Website: www.safetyzon.com