



## **Tilting Shaker TRTS-601**

[www.safetyzon.com](http://www.safetyzon.com) | [info@safetyzon.com](mailto:info@safetyzon.com)

## Description

Tilting Shaker TRTS-601 has a dual LCD screen that shows time and speeds separately. It has a working temperature range of 5 to 40 °C and a DC motor control system that controls the mechanism. It has a 7° inclination and a speed range of 10 to 80 revolutions per minute. A constant or countdown procedure is presented from 1 minute to 19h 59min. various inserts and fixings are included.

## Features:

- > It has a dual LCD screen, an LCD screen that displays time and speeds separately
- > The mechanism that is controlled by a DC motor control system
- > Provided with various inserts and fixings
- > Having a constant or countdown procedure mode

## Specifications :

<b>Angle of Inclination</b>	7°
<b>Data Access</b>	None
<b>Dimension (L×W×H)</b>	335×270×120 mm
<b>Maximum Load</b>	3kg
<b>Moisture</b>	0.8
<b>Motor Type</b>	DC motor
<b>Operation Model</b>	Time set, continuous operation
<b>Power</b>	20W
<b>Speed Indicator</b>	LCD
<b>Speed Range</b>	10 to 80rpm
<b>Time Range</b>	1 min to 19h59min
<b>Timer Indicator</b>	LCD
<b>Voltage</b>	100 to 240V, 50/60Hz
<b>Weight</b>	3.2kg
<b>Working Temperature</b>	5 to 40°



West 110th Street Overland Park KS 66210, USA  
Email: [info@safetyzon.com](mailto:info@safetyzon.com) | Website: [www.safetyzon.com](http://www.safetyzon.com)