



Water Distiller TRWD-603

www.safetyzon.com | info@safetyzon.com

Description

Water Distiller TRWD-603 is corrosion and aging resistant and has an intelligent automatic control. It is equipped with a stainless steel bowl with an electric heater. It condenses during the steaming process for the distilled water by electric heating. It has an output water level of 20 liters per hour.

Features:

- > It is corrosion-resistant and age resistant
- > Having a smart automatic control
- > Electric heating causes steam condensation into distilled water
- > It is outfitted with distillation generating equipment, which includes a condenser and other components
- > It consists of a stain-free steel electric double distiller

Specifications :

Output	20L/hour
Package	350×495×790 mm
Power	AC380V/13.5 KW
Weight	10 kg



West 110th Street Overland Park KS 66210, USA
Email: info@safetyzon.com | Website: www.safetyzon.com